



Microcone® Implant



Quattrocone Implant

NI	RI								
Ø 3,0 mm	Ø 3,5 mm	Ø 4,0 mm	Ø 4,5 mm	Ø 5,0 mm	Ø 3,5 mm	Ø 3,8 mm	Ø 4,3 mm	Ø 5,0 mm	
									
6,5 mm -	-	2-01-35	2-01-41	2-01-47	6,5 mm -	-	-	-	
7 mm -	-	-	-	-	7 mm -	-	3-01-16	3-01-06	-
8 mm -	2-01-30	2-01-36	2-01-42	2-01-48	8 mm -	-	-	-	-
9 mm -	2-01-31	2-01-37	2-01-43	2-01-49	9 mm 3-01-02	3-01-17	3-01-07	-	
11 mm 1-01-06	2-01-32	2-01-38	2-01-44	2-01-50	11 mm 3-01-03	3-01-18	3-01-08	3-01-13	
13 mm 1-01-07	2-01-33	2-01-39	2-01-45	2-01-51	13 mm 3-01-04	3-01-19	3-01-09	3-01-14	
15 mm 1-01-08	2-01-34	2-01-40	2-01-46	2-01-52	15 mm -	3-01-20	-	3-01-15	

Closure screw



Gingiva former



Auxiliary Instruments



Screw Driver short
0-13-22



Screw Driver long
0-13-23



Torque Ratchet 10-40 Ncm
0-13-28

Implant-level Impression

Implant pick-up
Open tray Impression

	NI	RI
short	-	2-04-01
long	1-04-01	2-04-02

Scanbody

	NI	RI
1-09-03		
2-09-10		

Temporary abutment

	NI	RI
1-17-04		2-17-07
1-06-01		2-06-03

Incl. screw

	NI	RI
1-05-01		2-05-01

Analogue

	NI	RI
1-05-01		2-05-01

CAD/CAM Analogue

	NI	RI
1-05-02		2-05-02

Implant-level Cementable/Screw retained Abutments

Titanium base
Including Screw
(2-06-03)

	RI	RI	RI	RI	RI
Chimney height	4,0 - 8,0 mm				
Gingiva height	1,0 mm	1,5 mm	2,0 mm	2,5 mm	3,0 mm
SSC Flex	2-09-26	-	2-09-27	-	2-09-28
SSC Flex Molar	-	2-09-29	-	2-09-30	-
SSC Flex Rotating	2-09-32	-	2-09-33	-	2-09-34
SSC Flex Rotating Molar	-	2-09-35	-	2-09-36	-

Recommended torque 25 Ncm

Titanium base
Including Screw
(1-06-01)

	NI	NI
Chimney height	4,0 - 8,0 mm	
Gingiva height	1,0 mm	2,0 mm
SSC Flex	1-09-12	1-09-13
SSC Flex Rotating	1-09-14	1-09-15

Recommended torque 15 Ncm

Abutment level Cementable/Screw retained Abutments

Multi-unit abutment

	RI	RI	RI		RI	RI	RI	RI	RI	RI
	0°				17°			30°		
Gingiva height	1,5 mm	3,5 mm	5,5 mm	-	2,5 mm	3,5 mm	5,5 mm	3,0 mm	4,0 mm	5,5 mm
Article No.	2-31-01	2-31-02	2-31-03	Type 1 Type 2	2-31-04	2-31-05	2-31-06	2-31-07	2-31-08	2-31-09

Recommended torque 25 Ncm

Multi-unit prosthetic components

	Titanium Cap	Titanium Base	Modelling sleeve
Article No.	0-31-06	0-31-09	0-31-11
Incl. Screw	0-31-02		

Recommended torque 15 Ncm

Multi-unit prosthetic accessories

	Analog	Implant pick up	Multi-unit cover cap	Scanbody	Analog CAD/CAM
Article No.	0-31-05	0-31-04	0-31-03	0-31-01	0-31-10

Placement instrument
Multi-unit abutment straight

0-13-76

MedentiBASE abutments

	RI	RI	RI	RI	RI
	0°				
Gingiva height	0,5 mm	1,5 mm	2,5 mm	3,5 mm	4,5 mm
Article No.	2-28-01	2-28-02	2-28-03	2-28-04	2-28-05

Recommended torque 25 Ncm

MedentiBASE prosthetic components

	Titanium Cap	Titanium Base	Plastic Cap
Article No.	0-28-06	0-28-05	0-28-08
Incl. Screw	0-06-03		

Recommended torque 15 Ncm

MedentiBASE accessories

	Analog	Implant pick up	Cover Cap	Scanbody	Analog CAD/CAM
Article No.	0-28-20	0-28-12	0-28-04	0-28-13	0-28-21

Placement instrument
MedentiBASE abutments

0-13-37

Novaloc®

Novaloc® Processing package titanium

2010.601

Novaloc® abutment

	RI	RI	RI		
	0°			15°	
Gingiva height		Type 1	Type 2		
1,0 mm	2-23-01	2-23-06	2-23-11		
2,0 mm	2-23-02	2-23-07	2-23-12		
3,0 mm	2-23-03	2-23-08	2-23-13		
4,0 mm	2-23-04	2-23-09	2-23-14		
5,0 mm	2-23-05	2-23-10	2-23-15		

Type 1 Type 2

Placement instrument
Novaloc® abutment angled

0-13-59