

Straumann®  
TLX Implant System


System Overview

# TLX SURGERY

∅ 3.75 mm 

∅ 4.5 mm 

∅ 5.5 mm 

∅ 6.5 mm 

Standard Plus (SP)  
RT ∅4.8mm




Standard Plus (SP)  
RT ∅4.8mm



Standard Plus (SP)  
WT ∅6.5mm



Standard Plus (SP)  
WT ∅6.5mm



Roxolid® SLActive®

6 mm	035.3106S
8 mm	035.3108S
10 mm	035.3110S
12 mm	035.3112S
14 mm	035.3114S
16 mm	035.3116S
18 mm	035.3118S

6 mm	035.3506S
8 mm	035.3508S
10 mm	035.3510S
12 mm	035.3512S
14 mm	035.3514S
16 mm	035.3516S
18 mm	035.3518S

6 mm	035.3706S
8 mm	035.3708S
10 mm	035.3710S
12 mm	035.3712S
14 mm	-
16 mm	-
18 mm	-

6 mm	035.3806S
8 mm	035.3808S
10 mm	035.3810S
12 mm	035.3812S
14 mm	-
16 mm	-
18 mm	-

Roxolid® SLA®

6 mm	035.2106S
8 mm	035.2108S
10 mm	035.2110S
12 mm	035.2112S

6 mm	035.2506S
8 mm	035.2508S
10 mm	035.2510S
12 mm	035.2512S

6 mm	035.2706S
8 mm	035.2708S
10 mm	035.2710S
12 mm	035.2712S

6 mm	035.2806S
8 mm	035.2808S
10 mm	035.2810S
12 mm	035.2812S

Closure Caps

RT

RT

WT

WT

RT Closure Cap


	∅ 2.7 mm	∅ 5.5 mm
 H 0 mm	036.3200S	-
H 1.5 mm	-	036.1201S

WT Closure Cap


	∅ 2.7 mm
 H 0 mm	036.3200S

Healing Caps




RT Healing Cap

	∅ 5.5 mm
 H 2.0 mm	036.1202S
H 3.0 mm	036.1203S
H 4.5 mm	036.1204S




WT Healing Cap

	∅ 7.2 mm
 H 2.0 mm	036.2202S
H 3.0 mm	036.2203S
H 4.5 mm	036.2204S

Important Information

	Prosthetic platform	Implant Diameters	Implant platform
RT ••	Two dots 	∅ 3.75 mm ∅ 4.5 mm	
			

Components laser-marked RT (Regular, TorcFit™) and with two dots can be used on all TLX implants with the RT neck.

	Prosthetic platform	Implant Diameters	Implant platform
WT •••	Three dots 	∅ 5.5 mm ∅ 6.5 mm	
			

Components laser-marked WT (Wide, TorcFit™) and with three dots can be used on all TLX implants with the WT neck.

Auxiliary Instruments



Torque Control Device  
066.1100



Ratchet  
046.119



Surgical Handle  
066.4000



RT/RN Bone Profiler 28 mm  
044.304




WT/WN Bone Profiler 25mm  
044.305



Guiding Pin for TLX Implant  
049.406S

RT - Regular TorcFit™  
Ø 4.8 mm



WT - Wide TorcFit™  
Ø 6.5 mm



Impression Post Open Tray	Impression Post Closed Tray	Implant Analog	Repositionable Implant Analog	CARES® Mono Scanbody
036.1000 036.1001	036.1002	036.0102	036.2100	036.3220

Impression Post Open Tray	Impression Post Closed Tray	Implant Analog	Repositionable Implant Analog	CARES® Mono Scanbody
036.2000	036.2002	036.1100	036.2102	036.3220

Temporary Abutment, crown	Temporary Abutment, bridge	Replacement Screw
037.1000	037.1001	036.3110

Temporary Abutment, crown	Temporary Abutment, bridge	Replacement Screw
037.2000	037.2001	036.3110

Variobase® for Crown	Variobase® for Crown AS	Variobase® for Bridge/Bar Cylindrical	Variobase® C
037.1201 037.1211	037.1203 037.1212	037.1204 037.1213	037.1205
Replacement Screw		Replacement Screw for Variobase® for Crown AS	
036.3110		036.3111	

Variobase® for Crown	Variobase® for Crown AS	Variobase® for Bridge/Bar Cylindrical	Variobase® C
037.2201 037.2211	037.2203 037.2212	037.2204 037.2213	037.2205
Replacement Screw		Replacement Screw for Variobase® for Crown AS	
036.3110		036.3111	

Novaloc Abutment							
	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm	2010.601-STM
0°	037.1500	037.1501	037.1502	037.1503	037.1504	037.1505	
15°	-	037.1510	037.1511	037.1512	037.1513	037.1514	

Novaloc Abutment							
	1 mm	2 mm	3 mm	4 mm	5 mm	6 mm	2010.601-STM
0°	037.2500	037.2501	037.2502	037.2503	037.2504	037.2505	
15°	-	037.2510	037.2511	037.2512	037.2513	037.2514	

Screwdriver 15mm 046.400	Screwdriver 21mm 046.401	Screwdriver 27mm 046.402

AS Screwdriver 15mm 046.786	AS Screwdriver 21mm 046.787	AS Screwdriver 27mm 046.788