

BLX IMPLANT

BLC IMPLANT

Closure Cap

Healing Abutment

AH = abutment height
GH = gingiva height
Ø = diameter



Auxiliary Instruments


Ø 3.5 mm 	Ø 3.75 mm 	Ø 4.0 mm 	Ø 4.5 mm 	Ø 5.0 mm 	Ø 5.5 mm 	Ø 6.5 mm 
RB (Regular Base) 	RB (Regular Base) 	RB (Regular Base) 	RB (Regular Base) 	WB (Wide Base) 	WB (Wide Base) 	WB (Wide Base) 
6 mm —	6 mm 061.4306	6 mm 061.5306	6 mm 061.6306	6 mm 061.7306	6 mm 061.8306	6 mm 061.9306
8 mm 061.3308	8 mm 061.4308	8 mm 061.5308	8 mm 061.6308	8 mm 061.7308	8 mm 061.8308	8 mm 061.9308
10 mm 061.3310	10 mm 061.4310	10 mm 061.5310	10 mm 061.6310	10 mm 061.7310	10 mm 061.8310	10 mm 061.9310
12 mm 061.3312	12 mm 061.4312	12 mm 061.5312	12 mm 061.6312	12 mm 061.7312	12 mm 061.8312	12 mm 061.9312
14 mm 061.3314	14 mm 061.4314	14 mm 061.5314	14 mm 061.6314	14 mm 061.7314	14 mm —	14 mm —
16 mm 061.3316	16 mm 061.4316	16 mm 061.5316	16 mm 061.6316	16 mm 061.7316	16 mm —	16 mm —
18 mm 061.3318	18 mm 061.4318	18 mm 061.5318	18 mm 061.6318	18 mm 061.7318	18 mm —	18 mm —

Ø 3.3 mm 	Ø 3.75 mm 		Ø 4.5 mm 		Ø 5.5 mm 	Ø 6.5 mm 
RB (Regular Base) 	RB (Regular Base) 		WB (Wide Base) 		WB (Wide Base) 	WB (Wide Base) 
6 mm —	6 mm 035.9206S		6 mm 035.9406S		6 mm 035.9706S	6 mm 035.9806S
8 mm 035.9008S	8 mm 035.9208S		8 mm 035.9408S		8 mm 035.9708S	8 mm 035.9808S
10 mm 035.9010S	10 mm 035.9210S		10 mm 035.9410S		10 mm 035.9710S	10 mm 035.9810S
12 mm 035.9012S	12 mm 035.9212S		12 mm 035.9412S		12 mm 035.9712S	12 mm 035.9812S
14 mm 035.9014S	14 mm 035.9214S		14 mm 035.9414S		14 mm 035.9714S	14 mm 035.9814S
16 mm 035.9016S	16 mm 035.9216S		16 mm 035.9416S		16 mm 035.9716S	16 mm 035.9816S
18 mm 035.9018S	18 mm 035.9218S		18 mm 035.9418S		18 mm —	18 mm —

RB (Regular Base)				WB (Wide Base)		
 H 0.4 mm 064.4100S				 H 0.5 mm 064.8102S		
				 No WB abutments on RB implants!		

Healing Abutments for Crown					Healing Abutments for Bridge				
AH 2 mm					AH 4 mm				
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm	GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	Ø 7.5 mm
1.5 mm				—	1.5 mm				—
2.5 mm					2.5 mm				
3.5 mm					3.5 mm				


Healing Abutments for Crown Wide Profile			Healing Abutments for Crown Wide Profile		
AH 2 mm			AH 4 mm		
GH	Ø 6.0 mm	Ø 7.5 mm	GH	Ø 6.0 mm	Ø 7.5 mm
0.75 mm			0.75 mm		
1.5 mm			1.5 mm		

 Torque Control Device 066.1100	 Ratchet 046.119	 BL Bone Profiler 1 026.0022	 BL Bone Profiler 2 026.0023	 BL Bone Profiler 3 026.0024	 Guiding Pin for BLX Implants 066.0025S	 RB/WB Abutment-Removal Screw 065.0007	 Removal Tool or BLX Basal Screw length 27 mm 065.0008	 Removal Tool for BLX Basal Screw length 21 mm 065.0009
--	---	---	--	---	--	---	---	--

Straumann® BLX Implant TorcFit™ Connection

RB/WB prosthetics

Fits on RB and WB implants (Ø 3.5 mm - Ø 6.5 mm)



WB prosthetics

Fits only on WB implants (Ø 5.0 mm - Ø 6.5 mm)



Impression Posts & Implant Analogs


Metal Scanbody & Repositionable Analog

Temporary Abutments


Impression Posts & Implant Analogs

Metal Scanbody & Repositionable Analog

Temporary Abutments



Variobase® XC & Burn-out Copings



Screw Retained Abutments & Novaloc

