



















**Straumann®
SLActive®
Roxolid®
BLX
IMPLANT**

Ø 3.3 mm 	Ø 3.75 mm 	Ø 4.0 mm 	Ø 4.5 mm 	Ø 5.0 mm 	Ø 5.5 mm 	Ø 6.5 mm 
RB (Regular Base) Ø 3.5 mm 	RB (Regular Base) Ø 3.75 mm 	RB (Regular Base) Ø 4.0 mm 	RB (Regular Base) Ø 4.5 mm 	WB (Wide Base) Ø 5.0 mm 	WB (Wide Base) Ø 5.5 mm 	WB (Wide Base) Ø 6.5 mm 
6 mm 061.3306	6 mm 061.4306	6 mm 061.5306	6 mm 061.6306	6 mm 061.7306	6 mm 061.8306	6 mm 061.9306
8 mm 061.3308	8 mm 061.4308	8 mm 061.5308	8 mm 061.6308	8 mm 061.7308	8 mm 061.8308	8 mm 061.9308
10 mm 061.3310	10 mm 061.4310	10 mm 061.5310	10 mm 061.6310	10 mm 061.7310	10 mm 061.8310	10 mm 061.9310
12 mm 061.3312	12 mm 061.4312	12 mm 061.5312	12 mm 061.6312	12 mm 061.7312	12 mm 061.8312	12 mm 061.9312
14 mm 061.3314	14 mm 061.4314	14 mm 061.5314	14 mm 061.6314	14 mm 061.7314	–	–
16 mm 061.3316	16 mm 061.4316	16 mm 061.5316	16 mm 061.6316	16 mm 061.7316	–	–
18 mm 061.3318	18 mm 061.4318	18 mm 061.5318	18 mm 061.6318	18 mm 061.7318	–	–

**Straumann®
SLActive®
Roxolid®
BLC
IMPLANT**

Ø 3.3 mm 	Ø 3.75 mm 	Ø 4.5 mm 	Ø 5.5 mm 	Ø 6.5 mm 
RB (Regular Base) Ø 3.5 mm 	RB (Regular Base) Ø 3.75 mm 	WB (Wide Base) Ø 4.5 mm 	WB (Wide Base) Ø 5.5 mm 	WB (Wide Base) Ø 6.5 mm 
6 mm 035.90085	6 mm 035.92065	6 mm 035.94065	6 mm 035.97065	6 mm 035.98065
8 mm 035.90105	8 mm 035.92085	8 mm 035.94085	8 mm 035.97085	8 mm 035.98085
10 mm 035.90105	10 mm 035.92105	10 mm 035.94105	10 mm 035.97105	10 mm 035.98105
12 mm 035.90125	12 mm 035.92125	12 mm 035.94125	12 mm 035.97125	12 mm 035.98125
14 mm 035.90145	14 mm 035.92145	14 mm 035.94145	–	–
16 mm 035.90165	16 mm 035.92165	16 mm 035.94165	–	–
18 mm 035.90185	18 mm 035.92185	18 mm 035.94185	–	–

















Closure Cap

 H 0.4 mm 064.41005	 H 0.5 mm 064.81025	  No WB abutments on RB implants!
--	--	---

Healing Abutment

AH = abutment height
GH = gingiva height
Ø = diameter



Healing Abutments for Crown				Healing Abutments for Crown				Healing Abutments for Bridge			Healing Abutments for Crown Wide Profile			Healing Abutments for Crown Wide Profile									
AH 2 mm				AH 4 mm				AH 2 mm			AH 4 mm			AH 2 mm			AH 4 mm						
GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	GH	Ø 4.0 mm	Ø 5.0 mm	Ø 6.5 mm	GH	Ø 5.0 mm	Ø 5.0 mm	GH	Ø 6.0 mm	Ø 7.5 mm	GH	Ø 6.0 mm	Ø 7.5 mm	GH	Ø 6.0 mm	Ø 7.5 mm				
1.5 mm				1.5 mm				1.5 mm			0.75 mm			0.75 mm			1.5 mm			1.5 mm			
2.5 mm	064.42025	064.42125	064.42225	2.5 mm	064.42035	064.42135	064.42235	2.5 mm	064.42055	064.42155	064.42255	1.5 mm	064.82015	064.82035	0.75 mm	064.82025	064.82045	1.5 mm	064.82125	064.82145	1.5 mm	064.82135	064.82155
3.5 mm	064.42065	064.42165	064.42265	3.5 mm	064.42075	064.42175	064.42275	3.5 mm	064.42325	064.42335	1.5 mm	064.82125	064.82145	1.5 mm	064.82135	064.82155							


Auxiliary Instruments

 Torque Control Device 066.1100	 Ratchet 046.119	 BL Bone Profiler 1 026.0022	 BL Bone Profiler 2 026.0023	 BL Bone Profiler 3 026.0024	 Guiding Pin for BLX Implants 066.00255	 RB/WB Abutment-Removal Screw 065.0007	 Removal Tool or BLX Basal Screw length 27 mm 065.0008	 Removal Tool for BLX Basal Screw length 21 mm 065.0009
--	---	---	--	---	--	---	---	--

Straumann® BLX Implant TorcFit™ Connection

RB/WB prosthetics

Fits on RB and WB implants (∅ 3.5 mm - ∅ 6.5 mm)



WB prosthetics

Fits only on WB implants (∅ 5.0 mm - ∅ 6.5 mm)



Impression Posts & Implant Analogs


Metal Scanbody & Repositionable Analog

Temporary Abutments

Impression Posts & Implant Analogs

Metal Scanbody & Repositionable Analog

Temporary Abutments



For Crown, AH 11 mm

	∅ 3.8 mm	∅ 4.5 mm	∅ 6.0 mm
GH 1.5 mm	064.4362	064.4372	—
GH 2.5 mm	064.4363	064.4373	064.4382
GH 3.5 mm	064.4364	064.4374	064.4383

For Bridge, AH 5.5 mm

	∅ 4.5 mm
GH 1.5 mm	064.4352



Temporary Abutments for Crown, AH 11 mm

	WB
GH 0.75 mm	064.4391
GH 1.5 mm	064.4392

Variobase® & Burn-out copings

Variobase® for Crown, AH 5.5 mm

	RB/WB	
	∅ 3.8 mm	∅ 4.5 mm
GH 1.5 mm	062.4934	062.4944
GH 2.5 mm	062.4935	062.4945
GH 3.5 mm	062.4936	062.4946

Variobase® for Crown AS, AH 5.5 mm

	RB/WB
	∅ 4.5 mm
GH 1.5 mm	062.4972

Variobase® for Bridge/Bar AH 3.5 mm

	RB/WB
	∅ 4.5 mm
GH 1.5 mm	062.4961

Variobase® Basal Screw

	Basal Screw	Basal Screw AS
	065.0036	065.0037

Variobase® for Crown, AH 5.5 mm

	WB	
	∅ 5.5 mm	∅ 6.5 mm
GH 0.75 mm	062.4953	—
GH 1.5 mm	062.4954	—

Variobase® for Crown AS, AH 5.5 mm

	WB
	∅ 5.5 mm
GH 1.5 mm	062.4971

Screw Retained Abutments & Novaloc

Protective cap

	RB/WB		
	∅ 4.6 mm		
Angle	0°	17°	30°
0.75 mm	—	—	—
1.5 mm	062.47225	—	—
2.5 mm	062.47235	—	—
GH 3.5 mm	062.47245	062.47335	062.47435
GH 4.5 mm	062.47255	062.47345	062.47445
GH 5.5 mm	—	062.47355	062.47455

Protective cap

	H 5.1 mm
	024.4323
	024.4324

Impression posts

	∅ 4.6 mm	
Closed-tray	025.0014	
Open-tray	025.0012	

Analog for Screw-retained Abutments

	∅ 4.6 mm
	023.4756

Scanbody

	025.0081
--	----------

Novaloc®

	RB/WB	
	∅ 3.8 mm	
Angle	0°	15°
1.5 mm	062.4501	—
2.5 mm	062.4502	062.4507
3.5 mm	062.4503	062.4508
GH 4.5 mm	062.4504	062.4509
GH 5.5 mm	062.4505	062.4510
GH 6.5 mm	062.4506	062.4511
GH 7.5 mm	—	062.4512

Processing Package titanium

	2010.601-NOV
--	--------------

Titanium Copings

	∅ 4.6 mm	∅ 4.6 mm
Crown	023.4747	024.0023
Bridge	023.4751	024.0024

Variobase® for Bridge/Bar

	∅ 4.6 mm
	023.0028

Occlusal Screw

	023.4763
--	----------

Repositionable Analog

	025.0008
--	----------

Impression cap

	2010.722-NOV
--	--------------

Analog

	2010.721-NOV
--	--------------